



Name: | DOB: | ID: | 04/06/04 14:39:15

Patient Temperature (°C)


Rectal Temp



Rectal Rate



Skin



Scalp



Other



Radiant Warmth




Set Cap

Cooling (hh:mm)


00:00

Cap (°C)



8 Hours | 24 Hours | Full

PAUSE



PLAYBACK

Status: Standby

Touch **Start** to begin a cooling treatment.

Start | Pause | Marker | Patient | Tools | ?

Name: | DOB: | ID: | 04/06/04 14:39:15

Patient Temperature (°C)

Rectal Temp



Rectal Rate



Skin




Scalp



Other



Radiant Warmth



Set Cap

Cooling



Cap (°C)



Status



Touch **Start** to begin a cooling treatment.

8 Hours

24 Hours

R
E
C
T
A
L

C
A
P

Manage Files

Settings

Languages

Update

System
{TBD}

Graphs

Close Patient

Shutdown

Help Topics

Contact

About

Only one button can be in a pressed state at a time (i.e., only one menu open at a time). All menus shown open here only for comparison of menu contents.

Start

Pause

Marker

Patient

Tools

?