


Specification



Olympic Vac-Pac® Size 40-C

Intended Use

The Olympic Vac-Pac Size 40-C is a patient positioning support intended for use during medical imaging procedures, such as magnetic-resonance imaging (MRI) scans.

Sizes and Dimensions

Size ¹	Catalog No.	Dimensions ²		Shape
		Imperial	Metric	
40-C	51640-C	21 x 76 in.	53.3 x 193.1 cm	Large narrow rectangle

1. More than ten sizes of Vac-Pacs are available for use as patient positioning supports during physical therapy and surgical procedures. Contact Olympic Medical for product and ordering information.
2. Measurements (width x length) of compressed Vac-Pac. Depth is dependent upon positioning and bead compression.

Environmental

Temperature

- Use: 65–85°F 18–29°C
- Shipping/Storage: -40–158°F -40–70°C

Product Material

- **Outer skin:** 20-mil polyvinyl chloride (PVC) film
- **Valves:** White PVC
- **Inner beads:** Expanded polystyrene, 0.023–0.047 mm

Note: Vac-Pacs are latex-free.

Repair Parts

Patch Kit	301350-T509
Includes repair patch, plug, and 5-in-1 connector	
Patch Kit	301350-K509
Includes label applicator squeegee, repair patch, plug, and 5-in-1 connector.	
Patch Only	301350-509
5-in-1 Connector	100197-01
Replacement Valve Plug	300932-503
Vac-Pump , hand vacuum pump	51650

Instructions for Use

Refer to the *Olympic Vac-Pac Instruction Manual for Surgical Positioning* for information and precautions when using, cleaning, and repairing the Size 40-C.

Cleaning

After each use wipe Vac-Pacs clean using a soap-and-water solution or a common liquid disinfectant (e.g., Cavicide®, Sporidicin®).

Service

Punctured Vac-Pacs may require patching or may be returned to factory for repair or trade-in discounts.

Warranty

Unconditional one-year warranty against all manufacturing defects.

CE Conformance

Vac-Pacs conform to the European Device Directives (EDD) for Class I products.

Product Disposal

Disposal of the Vac-Pac does not require any special precautions. Follow local waste disposal regulations.



OLYMPIC MEDICAL

5900 First Ave S · Seattle, WA 98108 USA
Toll-free: 1-800-426-0353 (US/Canada)
Phone: 206-767-3500 · Fax: 206-762-4200
Web: www.OlympicMedical.com

© 2005 Olympic Medical Corporation. All rights reserved.

Olympic Medical reserves the right to make changes without notice in design, specifications, and models. The only warranty Olympic Medical makes is the express written warranty extended on the sale or rental of its products.

51640-SPC-EN0524P