


Specification



Olympic Vac-Pac®


Intended Use

Olympic Vac-Pac products are patient positioning supports intended for use during physical therapy and surgical procedures.

Sizes and Dimensions

Size	Catalog No.	Dimensions ¹		Shape
		Imperial	Metric	
6	51606	30 x 36.5 in.	76.2 x 92.7 cm	Three sections / flaps
8	51608	34 x 52.5 in.	86.4 x 133.4 cm	Three sections / flaps
11	51611	16 x 12 in.	40.6 x 30.5 cm	U-shape
20	51620	20.8 x 11 in.	52.8 x 27.9 cm	Small rectangle
23	51623	18 x 20 in.	45.7 x 50.8 cm	Small U-shape cutout on 18-in. side, small
25	51625	36.5 x 14.8 / 18.5 in.	92.7 x 37.6 / 47.0 cm	Tapered rectangle
30	51630	28.7 x 35.5 in.	72.9 x 90.2 cm	Small U-shape cutout on 28.7-in. side, medium
31	51631	38.7 x 35.5 in.	98.3 x 90.2 cm	Small U-shape cutout on 38.7-in. side, large
32	51632	46.7 x 35.5 in.	118.6 x 90.2 cm	Small U-shape cutout on 46.7-in. side, extra large
35	51635	34.5 x 28 in.	87.6 x 71.1 cm	T-shape / flaps
38	51638	38.7 x 35.5 in.	98.3 x 90.2 cm	Shoulder positioner with U-shape and moon-shape cutouts
40	51640	30 x 76 in.	76.2 x 193.1 cm	Large rectangle

1. Measurements (width x length) of compressed Vac-Pac. Depth is dependent upon positioning and bead compression.



Environmental

Temperature

- Use: 65–85°F 18–29°C
- Shipping/Storage: -40–158°F -40–70°C

Product Material

- **Outer skin:** 20-mil polyvinyl chloride (PVC) film
 - **Valve:** White PVC
 - **Inner beads:** Expanded polystyrene, 0.023–0.047 mm
- Note:** Vac-Pacs are latex-free.

Repair Parts

- Patch Kit** 301350-T509
Includes repair patch, plug, and 5-in-1 connector
- Patch Kit** 301350-K509
Includes label applicator squeegee, repair patch, plug, and 5-in-1 connector.
- Patch Only** 301350-509
- 5-in-1 Connector** 100197-01
- Replacement Valve Plug** 300932-503
- Vac-Pump** (hand vacuum pump) 51650

Cleaning

After each use wipe Vac-Pacs clean using a soap-and-water solution or a common liquid disinfectant (e.g., Cavicide®, Sporidic®).

Service

Punctured Vac-Pacs may require patching or may be returned to factory for repair or trade-in discounts.

Warranty

Unconditional one-year warranty against all manufacturing defects.

CE Conformance

Vac-Pacs conform to the European Device Directives (EDD) for Class I products.

Product Disposal

Disposal of the Vac-Pac does not require any special precautions. Follow local waste disposal regulations.

